Kansai Medical University OMERO Server Construction

October 18, 2024 Ver. 1.1 Yasuhisa Motoyama, University Information Center Yousuke Nakano, Center for Clinical Anatomy

Build an OMERO server on Mac Studio ($\underline{M1 \text{ processor}}$). Describe the basic information and installation method in the construction.

♦Outline

In building the OMERO server on Mac Studio, there were few descriptions about installing it on MacOS, In addition, since <u>Mac Studio uses the M1 processor</u>, it was difficult to build the server using virtualization software such as UTM or VirtualBox due to architectural issues, so Docker was used to build the system.

Reference : OMERO Server Documentation

https://hub.docker.com/r/openmicroscopy/omero-web-standalone https://hub.docker.com/r/openmicroscopy/omero-server https://github.com/ome/omero-web-docker/issues/27

Prior information

· Contact the facility information officer for the following information.

IP address :

Subnet mask :

Default gateway : DNS1 : DNS2 :

Docker configuration (approximate)

Storage	:	1 TB (max. 2 TB)
Memory	:	112 GB (max. 128 GB)
CPU	:	20 cores (max. 40 cores)

Table of Contents

◇Pre	eparation for introduction	. 2
0.	Environmental Preparation (Docker installation)	. 2
⇔Co	nstruction of OMERO-Server	. 3
1.	Creating a container for a postgres database	. 3
2.	Obtain an image of OMERO- SERVER	. 3
3.	Creation of OMERO-SERVER container	. 3
⇔Co	nstruction of OMERO-Web	.4
1.	Obtain an image of OMERO-WEB	.4
2.	Creation of OMERO-WEB container	.4
⇔Co	nfirmation of OMERO start-up	. 5
1.	Confirmation of OMERO.insight activation	. 5
2.	Confirmation of OMERO.web activation	. 7

◇Preparation for introduction

0. Environmental Preparation (Docker installation)Download the application from the official Docker website.

©Docker installation
https://www.docker.com/products/docker-desktop/

Launch the <u>Docker</u> application. Account linkage is required, so use any email address to link accounts.

Launch a <u>terminal</u> and start <u>Docker</u>. ©command docker login

The "docker" command becomes available.

♦ Construction of OMERO-Server

- 1. Creating a container for a postgres database
- Creating a <u>postgres</u> container
 ©command
 docker run -d --name postgres -e POSTGRES_PASSWORD=postgres postgres
- 2. Obtain an image of OMERO- SERVER
 - \bigcirc command

docker pull --platform linux/amd64 openmicroscopy/omero-server

3. Creation of OMERO-SERVER container

 \bigcirc command (XXXX is set arbitrarily)

docker run -d --platform linux/amd64 --name omero-server --link postgres:db ¥

When ">" appears, enter the following

- -e CONFIG_omero_db_user=postgres ¥
- -e CONFIG_omero_db_pass=postgres ¥
- -e CONFIG_omero_db_name=postgres ¥
- -e ROOTPASS=<mark>XXXX</mark> ¥
- -р 4063:4063 -р 4064:4064 ¥

openmicroscopy/omero-server

%important%

M1 processor (ARM64 architecture), so the platform must be specified as amd64 with "<u>--platform</u> linux/amd64".

♦ Construction of OMERO-Web

- 1. Obtain an image of OMERO-WEB
 - \bigcirc command

docker pull openmicroscopy/omero-web-standalone --platform linux/amd64

Creation of <u>OMERO-WEB</u> container
©command (XXXX is set arbitrarily) docker run -d --name omero-web --platform linux/amd64 ¥ When ">" appears, enter the following --link omero-server:omero ¥ -e OMEROHOST=XXXX ¥ -p 4080:4080 ¥ openmicroscopy/omero-web-standalone

Confirmation of Container

©command docker ps -a --no-trunc

Ensure that three containers are created.



◇Confirmation of OMERO start-up

- 1. Confirmation of OMERO.insight activation
- Launch OMERO.insight

••••••••••••••••••••••••••••••••••••••				
localhost [Low]				
Username: root				
Password:				
Login Quit				
583 OMERO is distributed under the terms of the GNU GPL. For more information, visit openmicroscopy.org				
••• ••• OME				

Expand the setting screen from $^{\bigcirc}$ and click, Specify the IP address of the host terminal from $\textcircled{\bullet}$.

😧 Servers	\times				
Servers Enter a new server or select an existing one.	Z				
Server Address Examples my. server. name (default case) my. server. name:port (for non-default ports) 127. 0. 0. 1:1234 (if you have an IP address) wss://my. server. name (connecting via websocket) Connecting to					
localhost					
Add new server X					
Enter server address:					
172XXXXXXX					

 $\boldsymbol{\cdot}$ Log in with the root account specified when creating the container.

••••••••••••••••••••••••••••••••••••••
172.2 23.2.10 [Low]
Username: root
Password: ·····
Login Quit
5E3 OMEDO is distributed under the terms of the GNU GDL For more information, visit operativescopy ore
•••• OME

\cdot If you can log in with OMERO client, it's OK.

😧 root root connected to 172		
File Edit Window Help		
📧 💇 📰 🔳 😩 🍇 Display	/ Groups 💌	Search
Projects		
Attachments		
🔁 Tags		
Images	⊗	
Administration	S	

- 2. Confirmation of <u>OMERO.web</u> activation
- Access to OMERO.web

Start a web browser and enter "http://XXX.XXX.XXX.XXX:4080" in the URL to access the site.

"XXX.XXX.XXX.XXX.XXX" is the IP address of the host terminal.

If the following page is displayed, success.

	omero:4064 🗸
	Username:
	Password:
	Login
© 2 OMERO is distribute	OMERO.web 5.23.0. 2007-2022 University of Dundee & Open Microscopy Environment d under the terms of the GNU GPL. For more information, visit openmicroscopy.org Download OMERO.insight for Mac OS X, Windows, Linux

Confirm that you can log in.

		~		Search:		
👤 system root root 👻			General	Acquisition Preview		
Explore Tags Shares						
- 🔍 root root						
Orphaned Images						

Completed